



# Vac Pro

## Portable Suction Unit

The Vac series pump is APEX Medical's most recent generation of portable suction units and has been designed to combine function with aesthetic appeal and meet the needs of today's hygiene and infection control requirements. The Vac Pro is durable, lightweight and versatile in the provision of general suction usage.

# Vac Pro

## Portable Suction Unit

### Features ▼

- Smooth hygienic finish, water resistant, easy to clean and maintain.
  - Lightweight 2.8kg only.
  - Optional accessories available
- Single-patient use canister available in 800ml or 1800ml.  
Filter and Tubing.



Integrated Canister Holder



Portable design with ergonomic easy-to-carry handle.



Pressure regulator and easy-to-read analogue pressure gauge.



Detachable power cord allows you to quickly change to another cord to match local electrical outlet requirement.



Single-patient use canister available in 800ml or 1800ml.

### Specification

Electrical Rating	AC 110 - 120V, 60Hz
	AC 220 - 240V, 50 / 60Hz
Power Consumption	90W
Maximum Pressure	600 mmHg (adjustable)
Maximum Flow Rate	24 LPM
Net Weight	2.8 kg
Dimension	352(L) x 206(W) x 192(H) mm
Sound Level	< 55 dB

Distributed by

**Apex Medical Corp.**  
No. 9, Min Sheng St., Tu-Cheng,  
New Taipei City, 23679, Taiwan  
T +886 2 2268 3100  
F +886 2 2268 6525  
marketing@apexmedicalcorp.com  
www.apexmedicalcorp.com  
Print-2017/01 All rights reserved  
Specifications may change without notice.

